

PneuBack Chair™

SPINAL STRENGTHENING



The PneuBack Chair retrains muscles. Unweighting in certain postural altitudes with a fixed or stabilized pelvis allows for more effective lumbar distraction technique. This encourages the erector muscles to loosen up and allow more spinal mobilization.

STABILIZED UNWEIGHTING SYSTEM

- Train proprio and neuro receptors to begin recognizing and maintaining a new postural alignment at specific levels.
- Allows for strength evaluation at specific levels along the vertebral column.
- Isolates and strengthens specific levels along the vertebral column.
- Establishes specific strength and range of motion deficiencies.

Key Benefits

- Increase range of motion
- Stabilize the pelvis
- Postural changes
- Limited ROM
- Strength/recruitment issues
- Hyper/hypo mobile
- Post acute/chronic

Distributed by:

Pneumex
2605 North Boyer Ave
Sandpoint, Idaho 83864
800.447.5792



The PneuBack Chair with or without unweighting allows muscles to be rehabilitated.

"Every time my back went out I would spend days in bed. I tried everything short of surgery. The clinic in Sandpoint introduced me to the PneuBack Program. I am doing things I thought I would never do again. I ran with my dog for the first time in 15 years." Wes

"We love the PneuBack Program! We see a lot of back injury patients and most of them have already been to other clinics over their years of back pain. Once started in the program, patients often wish they could have come years ago." Don Gale, PT

Chair Studies on the Web

www.pneumex.com

Genuine Pneumex Equipment

www.pneumex.com

pneumex@pneumex.com

800.447.5792

Specifications

PneuBack Chair™

STRENGTHENING SYSTEM

ACCESSORIES INCLUDED:

PneuAngle™ Adjuster

Ability to un-weight in any vertical direction without moving Chair or PneuWeight station. Attaches to a PneuWeight station.



Knee Bolster

Used to keep legs in proper position while exercising in Chair



Chair Seat Extension

Provides extra height for shorter patients while exercising in Chair.



One-Cervical Halter

Specific stabilized mobilization while unweighted



Half Lumbar Roll

Provides support while exercising in Chair



Large Resistance Pad


Allows for more comfort during strengthening



CONFIGURATIONS AVAILABLE:



The PneuMAP generated protocols allow for precise exercise of spinal muscles.



Specifications

Height	Width	Length	Weight
63.25"	35.5"	66"	180lb

All PneuWeight systems require a compressor to operate.

Distributed by:

Pneumex
2605 North Boyer Ave
Sandpoint, Idaho 83864
800.447.5792



Genuine Pneumex Equipment

www.pneumex.com

pneumex@pneumex.com

800.447.5792

Rev 9-08-2014