

PneuVest™

UNWEIGHTING HARNESS



Medical Vectors Style

Our patented PneuVest is designed to comfortably support the patient during unweighting. Each vest also includes an buttock strap, leg straps and cuffs designed to alter and improve hip and gait problems.

Full Size Range

- Support for full range of unweighted activities.
- Easy to get in and out of, comfortable to wear.
- Unisex design adjusts to fit all body types.
- Design allows for normal diaphragm & rib interaction.
- Strapping options allow for gait correction.
- Made in U.S.A.



Back View With Straps-Cuffs



Left SI



Right SI



Enhance



Flatten

PneuVest	Description	Waist Size	Color
CHILD	XXSmall	<16" to 22"	Lt Blue
YOUTH	XSmall	22" to 28"	Berry
SMALL	Small	28" to 34"	Red
MEDIUM	Medium	34" to 44"	Blue
LARGE	Large	44" to 55"	Yellow
XLARGE	XLarge	55" to >67"	Purple

All Pneumex Vests are tested to 700 pounds load limit.

PneuVest Vectors Style vest has multiple adjustment vector points to help improve patient posture, increase or decrease lumbar lordosis, and unweight left or right SI independently.

Vectors has a sports bra type front and adjustable buckles on the back to assist with corrections.

Genuine Pneumex Equipment

www.pneumex.com

info@pneumex.com

800.447.5792