

# PneuWeight™

UNWEIGHTING SYSTEM



The PneuWeight Unweighting System allows patients to perform therapeutic activities in an upright, functional position without altering gait patterns in a safe and secure environment for both client and practitioner.

## A TRULY FLEXIBLE SYSTEM

- Unique cantilevered design for unobstructed access.
- The largest range of unweighting in the industry: **0 to 300 lbs.**
- 30" vertical oscillation — offering Balance and Exercise modes, as well as a drop-stop.
- 360 degree rotation allows multi-planar exercise.
- Patented vest with adjustability for posture and gait.
- Adjustable height range from 7'11" to 8'8".
- Sturdy and simple to use control panel.
- Customized pricing available.
- Available in single and double stations.
- Available with Analog control panel.



Patent 6672998

The system can be used over most exercise equipment, including but not limited to: treadmills, vibration platforms and balance equipment.

## Key Benefits

- Abilities patients
- Bariatric patients
- Orthopedic patients
- Strokes
- Back problems
- Gait training
- Cardiovascular
- Respiratory ailments
- Neurological patients
- Safety



### Analog Control Panel

Easy to read, easy to use analog control allows precise unweighting from 0 to 300 lbs. in 1 lb. increments. Features balance and exercise modes in familiar dial format.

All PneuWeight systems require a compressor to operate.

"I can't open my new clinic without the PneuWeight Unweighting System. It is one of the few capital purchases that I think is essential to my practice."

Dr. Bill O'Grady, DPT, MA, OCS,  
FAAOMPT, FAAPMa

"The PneuWeight system is the most effective program for dealing with treatment of back pain and prevention of injury. The results are consistent and reproducible between clinics."

Steve Sodorff, PT

Unweighting Studies on the Web

[www.pneumex.com](http://www.pneumex.com)

Genuine Pneumex Equipment

[www.pneumex.com](http://www.pneumex.com)

[pneumex@pneumex.com](mailto:pneumex@pneumex.com)

800.447.5792

# Specifications

# PneuWeight™

UNWEIGHTING SYSTEM

## ACCESSORIES

**Single Station = One PneuVest**  
**Double Station = Two PneuVest**

Light weight, easy to use.  
 Securely holds and lifts patients.

Four styles:

- Vectors
  - Sport
  - Abilities 3-in-1
  - X-Heart
- Sizes:  
 XSmall to  
 XLarge



Vectors  
Default PneuWeight



Sports  
Default TM-50



Abilities  
Default PneuWalker



X-Heart

### One-PneuGait Accessory Pack

1-Single Ankle Cuff



1-Double Ankle Cuff  
(for DropFoot)



1-PneuGait Foot  
Strap with Buckle



1-Gait Assist Band



### Your choice of style and size

#### One-Cable Assembly

Extra safety that allows  
 360° movement.



#### One-Chain Extender

Provides extra length on  
 cable Assembly



#### One-Safety Drop Stop

Limits the patient drop  
 to 2" if they start to fall.



## CONFIGURATIONS AVAILABLE:

## Specifications



Single Station

Includes One PneuVest



Double Station

Includes Two PneuVest

	Single	Double
Height	7'11" - 8' 8" 242 - 265 cm	7' 11" - 8' 8" 242 - 265 cm
Outside Base Width	4' 2" 127 cm	6' 4" 193 cm
Inside Base Width	3'8" 112 cm	3'8" 112 cm
Max Lifting	300 lb 136 k	300 lb* 136 k
* Each station		

All PneuWeight systems require a  
 compressor to operate.

### Distributed by:

Pneumex  
 2605 North Boyer Ave  
 Sandpoint, Idaho 83864  
 800.447.5792



Genuine Pneumex Equipment

www.pneumex.com

pneumex@pneumex.com

800.447.5792

Rev 8-8-2016